

mobi shower trolley



Safe showering

For many children and adults, conventional methods of bathing and showering can be difficult and even hazardous. The Mobi Shower Trolley is a powered, height-adjustable shower trolley specifically designed to enable the user to be safely and easily washed or showered by their carer.

Reduces back injury


The Mobi Shower Trolley offers a range of height adjustments to enable the carer to shower and change the user at a safe and suitable working height. This reduces the risk of back injury for the carer.

Stowaway side rails

The Mobi Shower Trolley has been designed with two unique safety side rails that can be easily stowed away beneath the trolley. This allows maximum clearance for safe, easy transfer and use.

mobi shower trolley

key features



- Mobile and easy to manoeuvre
- No installation required
- Stainless steel frame
- Drainage channel and pipe
- Soft comfortable mattress and wedge shaped pillow
- Available in 1 size
- Length of table is 1800mm
- Weight capacity of 220kg
- Travel range of 460mm
- Easy to clean padded surfaces
- Unique stowaway side rails
- Optional padded mattress as an accessory

Stainless Steel	Size	User Weight	Table Dimensions		Travel Range	Safety Sides		Safety Sides Operational Space	Castor Centres		Nominal Clearances	Shipping Weights
		kg	L mm	W mm		H mm	Width X mm		Height Y mm	Z mm		
3151	1	220	1800	750	465 - 925	840	195	200	1090	700	190	86

